





Ophthalmic Surgery

Ophthalmologie


NEEDLE-THREAD COMBINATIONS NADEL-FADEN-KOMBINATIONEN

Needle dimensions / Nadel-Abmessungen













		TYPE TYP	CIRCLE KREIS	CURVE KURVE	LENGTH LÄNGE	CHORD SEHNE	NEEDLE Ø NADEL-Ø
Cutting needle, lancetpoint 		HSL 7/15	1/2	175°	7.0 mm	5.0 mm	0.15 mm
		DSL 6/15	3/8	135°	6.2 mm	5.0 mm	0.15 mm
	Schneidende Nadel mit Lanzett-Spitze	DSL 6/20	3/8	135°	6.4 mm	5.0 mm	0.20 mm
Reverse cutting needle, micropoint 		DSM 7/20	3/8	110°	7.9 mm	6.4 mm	0.20 mm
		DSM 10/28	3/8	135°	11.0 mm	8.9 mm	0.28 mm
	Schneidende Nadel mit Mikro-Spitze	DSM 13/35	3/8	155°	13.0 mm	9.6 mm	0.35 mm
Scleral needle, sidecutting 		VSP 8/35	1/4	112°	7.9 mm	7.2 mm	0.35 mm
	Spatula Nadel, schneidend						

TYPE OF NEEDLE NADEL TYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	0.2	0.3	0.4	0.5	0.6
			LENGTH LÄNGE	10-0	9-0	8-0	7-0	6-0

 1/2 circle, lancetpoint,
cutting needle
1/2 Kreis, Lanzett-Spitze,
schneidende Nadel

2 HSL 7/15 	Nylon		30 cm	290 R 02				
---	-------	---	-------	----------	--	--	--	--

 3/8 circle, lancetpoint,
cutting needle
3/8 Kreis, Lanzett-Spitze,
schneidende Nadel

DSL 6/15 	Nylon		15 cm	74 R 02				
2 DSL 6/15 	Nylon		30 cm	274 R 02				
DSL 6/20 	Nylon		15 cm	95 R 02	95 R 03			
2 DSL 6/20 	Nylon		30 cm	295 R 02	295 R 03			
	Barraquer Silk/Seide		45 cm			295 H 04		

STERILE EO  0482


Packaging Unit: 12 Foils/Box
VPE: 12 Folien/Box







Ophthalmic Surgery

Ophthalmologie

NEEDLE-THREAD COMBINATIONS NADEL-FADEN-KOMBINATIONEN

TYPE OF NEEDLE NADEL TYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	0.4	0.5	0.7	1	1.5
			LENGTH LÄNGE	8-0	7-0	6-0	5-0	4-0

 3/8 circle, reverse cutting needle micropoint
3/8 Kreis, schneidende Nadel mit Mikrospitze

DSM 7/20 	Nylon	■	45 cm			81 R 07		
	Silk/Seide	■	45 cm		81 D 05	81 D 07		
2 DSM 7/20 	Silk/Seide	■	45 cm		281 D 05	281 D 07		
DSM 10/28 	Trofilene	■	45 cm			82 B 07		
	Nylon	■	45 cm			82 C 07		
	Silk/Seide	■	45 cm			82 D 07		
2 DSM 10/28 	Silk/Seide	■	45 cm			282 D 07		
DSM 13/35 	Trofilene	■	45 cm				83 B 1	83 B 15
	Silk/Seide	■	45 cm				83 D 1	
2 DSM 13/35 	Silk/Seide	■	45 cm					283 D 15

 1/4 circle, scleral side cutting needle
1/4 Kreis, Spatula-Nadel, schneidend

2 VSP 8/35 	Supramid	■	45 cm				298 U 1	298 U 15
	Polyester	■	45 cm				298 L 1	298 L 15

STERILE EO  0482

Packaging Unit: 12 Foils/Box
VPE: 12 Folien/Box