



BISTURI MONOUSO



MICROSURGERY KNIVES

STERILE · STAINLESS STEEL · SINGLE USE



Kai instruments are manufactured from state-of-the-art technologies accumulated through razor production and know-how based on blade production technologies for over 90 years.

- Polycarbonate plastic handle with color coding
- Transparent Polycarbonate plastic case
- Sterilized by E.O.G.
- Ready to use
- 5 pieces per dispenser box

Distinction of Kai's Stab knives
 Stab knife blades are made of razor blades with 3 bevel edge-setting. The razor edge-setting methods ensure the consistent sharpness and preciseness of stab knives.

Section of cutting edge

Stab knives

A-15F	Angle 15°	t=0.1mm	A-225F	Angle 22.5°	t=0.1mm
A-30F	Angle 30°	t=0.1mm	A-45F	Angle 45°	t=0.1mm

Crescent knives

CPU-26F	w=2.6mm	45° Bevel UP	CPD-26F	w=2.6mm	45° Bevel DOWN
CWU-26F	w=2.6mm	60° Bevel UP	CWD-26F	w=2.6mm	60° Bevel DOWN

Spoon knife

SPU-28F	w=2.8mm	45° Bevel UP

Paracentesis knife

PCT-10F	w=1.0mm



BISTURI MONOUSO

