



IMPIANTI RETINICI PER CHIRURGIA EPI-SCLERALE

Solid Silicone Implants

Bands, Laces and Sleeves

Circling Bands

Labtician offers a choice of three Circling Bands. All are designed to hold the accessory implant in place and to maintain proper pressure to form a permanent buckle without scleral erosion.

S 2969 (40 Style)

S 2987 (240 Style)

S 2950 (2950 Style)*

Silicone Lace™

The Silicone Lace makes Circling Band placement quicker and easier. The leader (which is removed at the end of the procedure) serves as a solid location to grasp the Band, and makes it easier to pass through a Silicone Sleeve.

S 5020 (40 Style)*

S 5025 (240 Style)*

Sleeves

Sleeves used to secure Circling Bands and Silicone Lace

S 3083 (3083 Style)*

S 3018 (70 Style)

S 3019 (270 Style)

Strips and Sleeves

Silicone Strips

Silicone Strips can be used alone, with wide grooved tires, or when trap door procedures are performed.

S 2970 (41 Style)

S 2971 (42 Style)

S 4050 (4050 Style)*

Sleeves

Sleeves used to secure Silicone Strips

S 3084 (3084 Style)*

S 3071 (72 Style)

Tires

Tire-shaped implants are designed for use under Circling Bands for wide scleral beds, for breaks near the ora serrata, for multiple breaks, and for high scleral buckles.

Convex

S 3010 (286 Style)

S 3014 (287 Style)

S 3014 L (287 WG Style)

S 3017 (288 Style)

S 3016 (289 Style)

Asymmetrical

S 2994 (275 Style)

S 2999 (276 Style)

S 2995 (278 Style)

S 2996 (280 Style)

S 2996 L (280 LG Style)

Concave

S 2986 (277 Style)

S 2997 (281 Style)

S 2991 (279 Style)

S 3004 (3004 Style)*

2.0 mm Concave

S 2977 (78 Style)

Wedges

Wedges are used to provide a greater radial buckling effect.

S 3021 (135 Style)

S 3022 (137 Style)

Grooved Strips

Grooved Strips are designed to be used under Circling Bands. Several designs are offered to fit the need for narrow or wide scleral beds, as well as high scleral buckling.

S 2965 (20 Style)

S 2967 (31 Style)

S 2968 (32 Style)

S 2992 (219 Style)

S 2998 (220 Style)

S 3013 (225 Style)

* Indicates products that are exclusive to Labtician.
 NOTE: Illustrations are not drawn to scale.



IMPIANTI RETINICI PER CHIRURGIA EPI-SCLERALE

Silicone Sponge Implants

Round	Part Thickness	Grooved	
S1982-1.5 (8215 Style) ★ 	S1984-2.5 (8425 Style) ★ 	S1983-2.3 (516G Style)	S1983-2.5 (510G Style) ★
S1982-2 (502 Style) 	S1984-3 (517 Style) 	S1983-2.8 (511G Style) ★ 	S1983-3 (506G Style)
S1982-2.5 (5025 Style) 	S1984-3.3 (515 Style) 	S1983-3.2 (519G Style) 	S1983-3.5 (509G Style)
S1982-3 (503 Style) 	S1984-5 (510 Style) 	S1983-4 (508G Style) 	S1983-5 (505G Style)
S1982-4 (504 Style) 	S1984-7.5 (511 Style) 	S1983-7.5 (507G Style) 	
S1982-5 (505 Style) 	Oblong	Oval	
S1982-7 (8270 Style) ★ 	S1986-3.5 (509 Style) 	S1981-4 (501 Style) 	S1981-5 (506 Style)
Scholda	S1986-4 (508 Style) 	S1981-7 (907 Style) 	S1981-7.5 (507 Style)
S1985-5 (8550 Style) ★ New 	Tunnel	Accessories	
Silicone Sponges are offered in many shapes which can be utilized for segmental buckling procedures.	S 1982-5T (505T Style) 	S 1987-7 (513 Style) 	S 1987-8 (514 Style)

