



## LENTI DIAGNOSTICHE



### OCULAR THREE MIRROR UNIVERSAL

Product Code	Image Mag	Static Fov	Lens Height	Contact Diameter
OG3M (Universal)	0.93x	140°	32mm	18mm
OG3M-2 (No-fluid)	0.93x	140°	32mm	16mm
OG3MF (With flange)	0.93x	140°	33mm	20mm
OG3MI (Universal)	0.93x	140°	28mm	15mm
OG3MP (Pediatric)	0.93x	140°	26mm	17mm
OG3MS (Small)	0.93x	140°	24mm	18mm
OG3MS-2 (Small no-fluid)	0.93x	140°	23mm	16mm
OG3M-13 (Small fissure no-fluid)	0.93x	140°	28mm	13mm



### OCULAR KARICKHOFF DIAGNOSTIC

Product Code	Image Mag	Static Fov	Lens Height	Contact Diameter
OJK	0.93x	140°	32.0mm	18mm
OJKF (No-fluid)	0.93x	140°	32.5mm	20mm



### OCULAR FUNDUS DIAGNOSTIC

Product Code	Image Mag	Static Fov	Lens Height	Contact Diameter
OGF	0.93x	36°	16.5mm	15.5mm
OGF-2 (No-fluid)	0.97x	35°	16.5mm	15.5mm



### OCULAR MAGNA VIEW GONIO

Product Code	Image Mag	Static Fov	Lens Height	Contact Diameter
OMVGL	1.3x	160°	23.5mm	15mm



### OCULAR SINGLE MIRROR GONIO

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diameter
OSMG	0.80x	170°	19.5mm	15mm



### OCULAR TWO MIRROR GONIO

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diameter
O2M	0.80x	170°	19.5mm	15mm
O2MG (No-fluid)	0.80x	170°	19.5mm	15mm



## LENTI DIAGNOSTICHE



### OCULAR THREE MIRROR 10MM GONIO

Product Code	Image Mag	Static Fov	Lens Height	Contact Diameter
OG3M-10	0.93x	140°	25mm	10mm



### OCULAR AUTOCLAVABLE THREE MIRROR

Product Code	Image Mag	Static Fov	Lens Height	Contact Diameter
OG3MAC-10	0.65x	150°	25mm	10mm



### OCULAR FOUR MIRROR MINI GONIO

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diameter	Ring Diameter
O4GF	0.80x	120°	22.5mm	15mm	23.5mm
O4GF-LR (Large ring)	0.80x	120°	26mm	15mm	32.5mm



### OCULAR MAXFIELD AC FOUR MIRROR GONIO

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diameter	Ring Diameter
O4MAC	0.61x	90°+	22mm	8.5mm	24.5mm
O4MAC-H (With handle)	0.61x	90°+	18mm	8.5mm	
O4MAC-LR (Large ring)	0.61x	90°+	28mm	8.5mm	31.5mm



### OCULAR GAASTERLAND FOUR MIRROR GONIO

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diameter	Ring Diameter
OG4MG	0.61x	90°+	22mm	8.5mm	24.5mm
OG4MG-H (With handle)	0.61x	90°+	18mm	8.5mm	
OG4MG-LR (Large ring)	0.61x	90°+	28mm	8.5mm	31.5mm



### OCULAR POSNER DIAGNOSTIC AND SURGICAL GONIOPRISM

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diam.	Handle Length
OPDSG (Round handle)	0.80x	80°	13mm	9mm	78.8mm
OPDSG-2 (Exagonal handle)	0.80x	80°	13mm	9mm	72.2mm
OPDSG-3 (Ergonomic handle)	0.80x	80°	13mm	9mm	92.8mm



## LENTI DIAGNOSTICHE



### OCULAR SUSSMAN FOUR MIRROR HAND HELD GONIOSCOPE

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diameter	Ring Diameter
OS4M	0.80x	80°	24.5mm	9mm	25mm
OS4M-2 (Large ring)	0.80x	80°	28.5mm	9mm	31.5mm



### OCULAR THORPE FOUR MIRROR GONIO

Product Code	Gonio Mag	Static Fov	Lens Height	Contact Diameter
OT4MG	0.80x	150°	32mm	18mm



### OCULAR KOEPPE DIAGNOSTIC

Product Code	Image Mag	Static Fov	Contact Diameter
OKL (Large)	1.50x	160°	19mm
OKM (Medium)	1.50x	160°	18mm
OKS (Small)	1.60x	160°	17mm