



LENTI PER YAG LASER



OCULAR ABRAHAM IRIDECTOMY

Product Code	Image Mag	Laser Spot Mag	Lens Height	Contact Diameter
OAIY	1.5x	0.67x	16.5mm	15.5mm



OCULAR ABRAHAM CAPSULOTOMY

Product Code	Image Mag	Laser Spot Mag	Lens Height	Contact Diameter
OAYA	1.8x	0.56x	16.5mm	15.5mm



OCULAR PEYMAN G. CAPSULOTOMY

Product Code	Image Mag	Laser Spot Mag	Lens Height	Contact Diameter
OPYG-12/12	1.8x	0.56x	16.5mm	15.5mm



OCULAR MANDELKORN IRIDOTOMY/CAPSULOTOMY

Product Code	Image Mag	Laser Spot Mag	Lens Height	Contact Diameter
OMIC	1.2x	0.83x	16.5mm	15.5mm



OCULAR MAGNA VIEW GONIO

Product Code	Gonio Mag	Laser Spot Mag	Static Fov	Lens Height	Contact Diameter
OMVGL	1.3x	0.77x	160°	23.5mm	14.5mm



OCULAR POLLACK IRIDOTOMY/GONIO

Product Code	Image Mag	Laser Spot Mag	Lens Height	Contact Diameter
OPIG	1.5x	0.65x	21mm	15mm



OCULAR PEYMAN WIDE FIELD

Product Code	Image Mag	Laser Spot Mag	Lens Height	Contact Diameter
OPY-12.5 (Anterior chamber)	1.40x	0.71x	16.5mm	15.5mm
OPY-18 (Mid vitreous)	1.41x	0.71x	16.5mm	15.5mm
OPY-25 (Deep vitreous)	1.36x	0.74x	14.7mm	16mm