



LENTI PER OFTALMOSCOPIA INDIRECTA



OCULAR MAXLIGHT HIGH MAG 78D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-HM	0.93x	1.07x	84°	139°	8mm	29.1mm



OCULAR MAXLIGHT STANDARD 90D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-STD	0.75x	1.34x	94°	153°	5mm	19.2mm
OI-STD-LR (Large ring)	0.75x	1.34x	94°	153°	5mm	19.2mm



OCULAR MAXLIGHT ULTRA MAG 60D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-UM	1.15x	0.87x	76°	131°	11mm	30mm



OCULAR MAXLIGHT 14D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Working Distance	Clear Aperture
OI-14	4.29x	0.23x	37°	72mm	52mm



OCULAR MAXLIGHT 20D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Working Distance	Clear Aperture
OI-20	2.97x	0.34x	50°	47mm	48mm



OCULAR MAXLIGHT TRIPLE TWO PANFUNDUS

Product Code	Image Mag	Laser Spot Mag	Static Fov	Working Distance	Clear Aperture
OI-222	2.72x	0.37x	60°	39mm	52mm



LENTI PER OFTALMOSCOPIA INDIRECTA



OCULAR MAXLIGHT 28D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Working Distance	Clear Aperture
OI-28	2.13x	0.47x	58°	29mm	38.2mm



OCULAR MAXFIELD 20D / OCULAR MAXAC 20D AUTOCLAVABLE

Product Code	Image Mag	Laser Spot Mag	Static Fov	Working Distance	Clear Aperture
OI-20M	2.97x	0.34x	50°	47mm	48mm
OI-20A (Autoclavable)	3.03x	0.33x	50°	47mm	48mm



OCULAR MAXFIELD 28D / OCULAR MAXAC 28D AUTOCLAVABLE

Product Code	Image Mag	Laser Spot Mag	Static Fov	Working Distance	Clear Aperture
OI-28M	2.11x	0.47x	58°	27mm	38.2mm
OI-28A (Autoclavable)	2.15x	0.47x	59°	28mm	38.2MM



OCULAR MAXFIELD 54D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-54M	1.10x	0.90x	86°	137°	9.8mm	29.1mm



OCULAR MAXFIELD 60D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-60M	0.99x	1.01x	88°	154°	9.8mm	29.1mm



OCULAR MAXFIELD 66D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-66M	0.91x	1.10x	91°	144°	7.8mm	27mm



LENTI PER OFTALMOSCOPIA INDIRECTA



OCULAR MAXFIELD 72D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-72M	0.83x	1.20x	102°	155°	6.8mm	27mm



OCULAR OSHER MAXFIELD 78D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-78M	0.77x	1.30x	98°	155°	7mm	27mm



OCULAR MAXFIELD 84D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-84M	0.71x	1.41x	105°	158°	5mm	27mm



OCULAR MAXFIELD 100D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-100M	0.60x	1.67x	110°	146°	4mm	29.1mm



OCULAR MAXFIELD 120D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-120M	0.50x	2.00x	120°	173°	4mm	21mm



OCULAR ULTRA VIEW SMALL PUPIL

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-SP	0.45x	2.22x	99°	152°	4mm	16mm



OCULAR MAXFIELD STANDARD 90D

Product Code	Image Mag	Laser Spot Mag	Static Fov	Dynamic Fov	Working Distance	Clear Aperture
OI-STDM	0.75x	1.34x	94°	153°	5mm	19.2mm
OI-STDM-LR (Large ring)	0.75x	1.34x	94°	153°	5mm	19.2mm