



LENTI VITRECTOMIA WIDE ANGLE



OCULAR LANDERS WIDE FIELD VITRECTOMY LENS

Product Code	Image Mag	Static Fov	Dynamic Fov	Lens Height
OLIV-WF	0.38x	130°	146°	12mm



OCULAR LANDERS EQUATORIAL II VITRECTOMY LENS

Product Code	Image Mag	Static Fov	Dynamic Fov	Lens Height
OLIV-EQ-2	0.65x	101°	131°	14.5mm



OCULAR WOLDOFF HIGH MAGNIFICATION VITRECTOMY LENS

Product Code	Image Mag	Static Fov	Dynamic Fov	Lens Height
OWIV-HM	.90x	57°	100°	13.5mm



OCULAR LANDERS NON-AUTOCLAVABLE WIDE FIELD VITRECTOMY LENS

Product Code	Image Mag	Static Fov	Dynamic Fov	Lens Height
OLIV-WFNA	0.38x	130°	146°	12mm



OCULAR LANDERS NON-AUTOCLAVABLE EQUATORIAL VITRECTOMY LENS

Product Code	Image Mag	Static Fov	Dynamic Fov	Lens Height
OLIV-EQNA	0.65x	101°	131°	14.5mm



OCULAR WOLDOFF NON-AUTOCLAVABLE HIGH MAGNIFICATION VITRECTOMY LENS

Product Code	Image Mag	Static Fov	Dynamic Fov	Lens Height
OWIV-HMNA	0.90x	57°	100°	13.5mm